

Table S2. Two cleaved amplified polymorphic sequences (CAPS) markers used to genotype Foz1 x Luc1 and Luc1 x Jer1 F₁ plants. Markers are adapted from Barbieri *et al.* 2012.

Cross	Marker	Primers (forward; reverse)	Temp. ¹	Restriction enzyme
Foz1 x Luc1	Bradi4g10930	GGTTGAGTTATTGCCGTCAG; CGCCTTTGGATCTACCAGTC	56°C	MseI
Luc1 x Jer1	Bradi4g10930	GGTTGAGTTATTGCCGTCAG; CGCCTTTGGATCTACCAGTC	56°C	MboI

¹Annealing temperature